


Allura Xper FD20 Floor

Preferred Room Layout

The layout shown below is based upon a typical equipment configuration with full movement range and should be considered as a general design guideline. Site conditions, application requirements, customer preferences, and/or equipment configuration may significantly impact suite design and equipment layout. It is recommended to request site-specific drawings from an Amber representative early in the design process.



Equipment Layout

Recommended Ceiling Height: 9' - 6³/₁₆" (2900mm)

Ceiling Heights other than recommended may impact equipment functionality, consult with Philips.

